


Medica Display for Medical Device

MEDDP-632F

32" Full-UHD Medical LCD Monitor



Features

- 32" Display with Full HD (1920 x 1080) native resolution
- Thin and compact design with impact resistant screen
- Fanless and ventless design, easy-to-clean
- Wide view angle 178°(H)/178°(V)
- True Flat Screen with PCT Touch
- Physical bottom OSD Design

Specifications

DISPLAY

Size	32"
Resolution	1920 x 1080 pixels
Max. Colors	16.7M Colors
Brightness (Typ.)	400 nits
Viewing Angle (CR=10)(Typ.)	178°(H)/178°(V)
Contrast Ratio (Typ.)	4000:1
Touch Screen	Projected Capacitive Multi Touch Screen

SYSTEM

OS Support	Windows® 7, Windows® 10
------------	-------------------------

I/O

Input Signal	VGA x 1, HDMI 1.4 x 1
USB Port	USB 2.0 x 1 for Touch Screen
Speaker	Built in 5W speaker x 2

MECHANICAL AND ENVIRONMENTAL

Operating Temperature	0°C to 35°C
Storage Temperature	-20°C to 60 °C
Storage Humidity	10% to 90% RH, noncondensing
VESA	100mm x 200mm and 200mm x 200mm
IP Rating	Front IP65
Power Requirement	External power adapter, AC IN: 100 ~ 240V, 12V
Dimension	799 x 478 x 63.7 mm
Package Size	980x260x700mm
Gross Weight	20 Kg
Net Weight	18 Kg
Certifications	CE: EN 60601-1-2:2015(V4.0), EN 60601-1:2006/A1:2013/A12:2014 (V3.1) FCC: Part 18 Class B

Ordering Information

- MEDDP-632F-P1-A1-0010
32" FHD.400nits LCD Monitor..HDMI x1, VGA x1.P-CAP.Med Adapter, Power cord