

Medical Power Panel PC

# 12"~ 17" Mini Mobile Computer Accessory Selection

## OPM-P05T

Venus-123/153/173 Battery Packs (2x)



## Specifications

Battery Type	253P / 6 Cell Lithium-ion battery
Battery Voltage	7.2V
Battery Capacity	4545 mAh / 32.72Wh
Warranty	1 year
Certification	UL,CE,IEC-60950-1,IEC-62133-1

### MECHANICAL AND ENVIRONMENTAL

Color	White
Dimension	127(L) x 60(w) x 21(H) mm
Net Packing Size	199(L) x 175(w) x 42(H) mm
Weight	305g/pcs
Gross Weight	0.9kg(1.98 lbs)
Operating Temperature	0°C ~ 40°C(32°F ~ 104°F)
Storage Temperature	-20°C ~ 60°C(-4°F ~ 140°F)

## Ordering Information

- OPM-P05T-A4  
Battery Kits (2x) 2S3P/7.2V for Venus-123/153/173

## UP-M62

Venus-123/153/173 6-slots Charger



## Optional Accessories

Part Number	UP-M62
Description	Battery Charger
Input Power	24V
Installations	Stationary
IO	1x Micro USB (Client) 1x DC-In
LED Indicator	Flickering Blue: Battery Charging Steady Blue: Battery Fully Charged Flickering Purple Green: Error
Dimensions	245.95 x 124.64 x 143.65mm
Weight	860g
Operating Temperature	-0°C ~ 40°C(32°F ~ 104°F)
Storage Temperature	-20°C ~ 60°C(-4°F ~ 140°F)

## Ordering Information

- UP-M62-A1-1010  
6-slots Charger for Venus-123/153/173