


Medica Display for Medical Device

MEDDP-615

15.6" Slim Medical Display



Features

- 15.6" 1920 x 1080 Resolution Display
- Projected Capacitive Multi-touch
- VGA input and HDMI input
- Stylish and elegant design

Specifications

DISPLAY

Size	15.6"
Resolution	Full HD 1920 x 1080
Max. Colors	16.7M
Brightness (Typ.)	300 nits
Viewing Angle (CR=10)(Typ.)	85°(H)/85°(V)
Contrast Ratio (Typ.)	700:1
Touch Screen	Projective capacitive touch screen

SYSTEM

OS Support	Windows® XP, Windows® 7, Windows® 10
------------	--------------------------------------

I/O

Audio	Audio-in x 1
USB Port	USB 2.0 Type B x 1 for Touch Screen
Signal	VGA x 1, HDMI x 1
Speaker	1 W speaker x 2 on back side

MECHANICAL AND ENVIRONMENTAL

Power Consumption	Full Loading :24.37Watts
VESA	75 x 75mm & 100 x 100mm
Power Requirement	DC12V / 5A power input
Operating Temperature	5°C ~ 40°C(41°F ~ 104°F)
Storage Temperature	-10°C ~ 60°C(14°F ~ 140°F)
Storage Humidity	10%~90%@35°C, non-condensing
Dimension	387.8 x 232.9 x 38.4 mm
Packing Size	425.3 x 269.2 x 45.6 mm
Gross Weight	6.3 kg
Net Weight	4.3 kg

Certifications	CE: EN 60601-1-2:2015(V4.0), EN 60950-1 EN 60601-1:2006/A1:2013/A12:2014 (V3.1) EN 62366:2007 + A1:2014 (V1.1) FCC: Part 18 Class B, Part 15 Class B UL: ANSI/AAMI ES60601-1:2012 (V3.1) cUL: CAN/CSA-C22.2 No. 60601-1:2014 (V3.1)
----------------	--

Ordering Information

- MEDDP-615HPN-A1-1010
15.6" LCD Monitor.300 nits.P-CAP touch,VGAx1.HDMIx1.DC 12V Adapter