




Medica Display for Medical Device

# MEDDP-722

21.5" Slim Medical Display



## Features

- 21.5" LCD Panel with LED Backlight
- Compact and TerseDesign
- DVI x 1, HDMI x 1, DPx1
- True Flat Screen with Capacitive Touch
- Speaker -2 x built-in
- Support 2M Camera

## Specifications

### DISPLAY

Size	21.5"
Resolution	FHD 1920 x 1080
Max. Colors	16.7M
Brightness (Typ.)	250 nits
Viewing Angle (CR=10)(Typ.)	178°(H)/178°(V)
Contrast Ratio (Typ.)	1000:1
Touch Screen	Projective capacitive touch screen (optional)
Power Requirement	External power adapter, AC IN: 100~240V, DC 12V/3A

### SYSTEM

OS Support	Windows® 10, Windows® 11
Webcamt	2MP

### I/O

Input Signal	DVI-D x1, HDMI 1.4 x1, VGA x1
USB Port	USB 2.0 Type B x1 for Touch Screen
Audio	Line-in x1
Speaker	Built in 5W speaker x2

### MECHANICAL AND ENVIRONMENTAL

Operating Temperature	5°C ~ 40°C(41°F ~ 104°F)
Storage Temperature	-10°C ~ 60°C(14°F ~ 140°F)
Storage Humidity	10%~90%@35°C, non-condensing
VESA	75mm x 75mm, 100mm x 100mm
Dimension (W/H/D)	546 x 351 x 56 mm
Packing Size (W/H/D)	644 x 489 x 211 mm
Gross Weight	8.4kg (18.52lb)
Net Weight	5.7kg (12.57lb)
Certifications	CE: EN 60601-1-2:2015 + A1:2021 (V4.1), EN 60601-1:2007/A1:2013/A2:2021 (V3.2) FCC: Part 15 Subpart B / Part 18 Subpart C UL: ANSI/AAMI ES60601-1: 2005 & A1:2012 & A2:2021 (V3.1) cUL: CAN/CSA-C22.2 No. 60601-1:2014 (V3.1) CCC

## Optional Accessories

- **125530036A**  
(AOH)(TF)Medical Power Adaptor.AC/DC.100 ~ 240V.12V.3A.36W. Adapter.ATM036T-P120(Level VI)
- **1700252001**  
(AOH)(TF)DVI-D Cable.(24+1)P DVI (M).(24+1)P DVI (M).200cm
- **1704191800**  
(AOH)(TF)HDMI Cable.HDMI 19P(F) Molding.HDMI 19P(M) Molding.180cm
- **1704151800**  
(AOH)(TF)VGA Cable.D-Sub 15P(M) Molding.D-Sub 15P(M) Molding.180cm
- **1700031802**  
(AOH)(TF)Aidio Cable.3.5mm 3Pole Stereo Plug(Black).3.5mm 3Pole Stereo Plug(Green).180cm

## Ordering Information

- **MEDDP-722-G1-A1-0010**  
Monitor. 22".250 nits.DC 12V.Glass
- **MEDDP-722-G1-A2-0010**  
Monitor. 22".250 nits.DC 12V.Glass.2M Camera
- **MEDDP-722-P1-A1-0010**  
Monitor. 22".250 nits.DC 12V.PCT
- **MEDDP-722-P2-A1-0010**  
Monitor. 22".250 nits.DC 12V.PCT.2M Camera