



Medica Display for Medical Device

MEDDP-822

21.5" Slim Medical Display



Optional Accessories

- 125530036A
(AOH)(TF)Medical Power Adaptor.AC/DC.100 ~ 240V.12V.3A.36W.Adapter. ATM036T-P120(Level VI)
- 1700252001
(AOH)(TF)DVI-D Cable.(24+1)P DVI (M).(24+1)P DVI (M).200cm
- 1704191800
(AOH)(TF)HDMI Cable.HDMI 19P(F) Molding.HDMI 19P(M) Molding.180cm
- 1704151800
(AOH)(TF)VGA Cable.D-Sub 15P(M) Molding.D-Sub 15P(M) Molding.180cm
- 1700031802
(AOH)(TF)Aidio Cable.3.5mm 3Pole Stereo Plug(Black).3.5mm 3Pole Stereo Plug(Green).180cm

Features

- 21.5" LCD Panel with LED Backlight
- High Brightness 500nits
- Compact and TerseDesign
- DVI x 1, HDMI x 1, VGAX1
- True Flat Screen with Capacitive Touch
- Speaker built-in x 2
- Physical OSD key bottom

Specifications

DISPLAY

Size	21.5"
Resolution	FHD 1920 x 1080
Max. Colors	16.7M
Brightness (Typ.)	500 nits
Viewing Angle (CR=10)(Typ.)	178°(H)/178°(V)
Contrast Ratio (Typ.)	1000:1
Touch Screen	Projected Capacitive Multi Touch Screen
Power Requirement	External power adapter, AC IN: 100~240V, DC 12V/3A

SYSTEM

OS Support	Windows® 10, Windows® 11
------------	--------------------------

I/O

Input Signal	DVI-D x1, HDMI 1.4 x1, VGA x1
USB Port	USB 2.0 Type B x1 for Touch Screen
Audio	Line-in x1
Speaker	Built in 5W speaker x 2

MECHANICAL AND ENVIRONMENTAL

Operating Temperature	0°C ~ 35°C(32°F ~ 95°F)
Storage Temperature	-20°C ~ 60°C(-4°F ~ 140°F)
Storage Humidity	10% to 90% RH, noncondensing
VESA	75mm x 75mm, 100mm x 100mm
Dimension (W/H/D)	526 x 324 x 54mm
Packing Size (W/H/D)	645 x 378 x 210 mm
Gross Weight	6.8kg (14.99lb)
Net Weight	5.0kg (11.02lb)
Certifications	CE: EN 60601-1-2:2015 + A1:2021 (V4.1), EN 60601-1:2007/A1:2013/A2:2021 (V3.2) FCC: Part 15 Subpart B / Part 18 Subpart C UL: ANSI/AAMI ES60601-1: 2005 & A1:2012 & A2:2021 (V3.1) cUL: CAN/CSA-C22.2 No. 60601-1:2014 (V3.1)

Ordering Information

- MEDDP-822-G1-A1-0010
Monitor. 22".500 nits.DC 12V.Glass.HDMI cable.Medical Adapter
- MEDDP-822-P1-A1-0010
Monitor. 22".500 nits.DC 12V.P-CAP.HDMI cable.Medical Adapter