



Medical Computer for Medical Device

MEDPC-9200

High Performance Medical Grade PC with Intel® 6th Generation Core™ iSeries



Features

- Intel® 6/7th Generation Core i, Processor
- Supports DDR4 up to 32GB
- Extension area: PCI-Express [x1] x 1(optional)
- Two Independent , Display: DP/HDMI
- Supports TPM 2.0
- Modularize Design
- Support Inte Modivus AI card

Application

- HIS, LIS, RIS
- Pharmaceutical industry
- Biotech Lab
- Signage in hospital
- Equipment control

Specifications

MAIN SPECIFICATIONS

Processor	Intel® 7th/6th Generation Core™ i7/i5/i3
System Memory	Supports Dual Channel DDR4 SODIMM up to 32 GB
OS Support	Windows® 7, Windows® 8, Windows® 10, Linux®
Extension area	PCI-Express [x1] x 1(optional)
Storage Disk Drive	2.5" SATA SSD x 1
Security	TPM 2.0
Power Requirement	DC 12V

I/O


USB	USB 3.0 x 4 , USB 2.0 x 2
Serial Port	RS-232 x 2
Ethernet	Gigabit LAN x 2
Video Out	Display Port x 1 (Support up to 1920 x 1200 @60 Hz) , HDMI 1.4 x 1
Audio	Mic-in, Line-out

MECHANICAL AND ENVIRONMENTAL

Operating Temperature	0°C ~ 40°C(32°F ~ 104°F)
Storage Temperature	-20°C ~ 60°C(-4°F ~ 140°F)
Dimension	197 x 197 x 59.5 (mm) for main layer 197 x 197 x 44 (mm) for extension layers
Net Weight	2.5 kg in the single layer, 3.7 kg in two layers
Package Size	300 x 280 x 95 (mm)
Gross Weight	3.3 kg (7.25 lb) in the single layer, 4.0 kg (7.25 lb) in two layer
Certifications	CE : EN 60601-1-2: 2015 + A1: 2021 (V4.1), EN 60601-1: 2007 + A1:2013 +A2:2021 (V3.2) EN 55032:2015/A1:2020 EN55035:2017/A11:2020(ITE) IEC62368-1:2020+A11:2020 (ITE) FCC : Part 15B/ Part 18 UL : ANSI AAMI ES60601-1:2005/A1:2012/A2:2021(V3.2) cUL : CAN/CSA-C22.2 No. 60601-1:14 (IEC 60601-1:2005+A1:2012+A2:2020,MOD) (V3.2)

High Performance Medical Grade PC with Intel® 6th Generation Core™ iSeries

Dimension / Unit: mm



MEDPC-9200



Ordering Information

- MEDPC-9200-N1-A1-0010
Medical Box PC.i7-6700TE.Fanless.DC 12V
- MEDPC-9200-N2-A1-0010
Medical Box PC.i5-6500TE.Fanless.DC 12V

Optional Accessories

- 1255300907
(AOH)(TF) Power Adapter.100-240Vac.50/60Hz.1.5A..12Vdc.7.5A.90W.FSP.FSP090M-RHA(w/Lock)
- OPM-C16W-A7
WLAN Kit.802.11a/b/g/n w/BT 4.0.Atheros AR5B22
- M92-SEC-A1
MEDPC-9200.extension area with PCIe[x1] riser
- OPM-T029-A0
Intel modivius AI card for MEDPC-9200
- OPM-C15W-A7
WLAN Kit.802.11a/b/g/n/ac BT 5.0 Intel AC9260