


Medical Power Bank

# OPM-P02T-00

Li-ion 3S3P Battrey Pack 10.9V 8400mAh



## Features

- Compatible with all Venus-series and all UPower-series
- X1.5 battery capacity increase
- Supports remote management software (ORION)
- Battery capacity indicator button
- Medical grade certified (with Venus)
- Lightweight design

## Specifications

### MAIN SPECIFICATIONS

Battery Type	3S3P / 9 Cell Li-ion battery
Battery Voltage	10.9V
Battery Capacity	8400mAh
Warranty	1 year
Charging time	3hrs for 90% ; 3.5hrs for 100% ( with UPower Pro-series)
Safety	Supports protection from over-voltage (input and output), over-current (input and output), short circuit, over-charge, over-discharge, and battery Positive Temperature Coefficient(PTC)
Compatible	Venus all series, UPower all series
Certification	UL, CE, IEC-60950-1, IEC-62133-1

### MECHANICAL AND ENVIRONMENTAL

Architecture	ABS + PC
Color	White
Dimension	122(L) x 79(W) x 50 (H) mm
Net Weight	0.5kg (1.10 lbs)
Packing Size	220(L) x 180(W) x 80(H) mm
Gross Weight	1.3kg (2.87 lbs)

## Ordering Information

- **9789BAXL14**  
2XL Battery Kit. Li-ion 3S3P, 10.9V, 8400mAh. Volume Order (via Ocean)
- **9789BAXL14**  
2XL Battery Kit. Li-ion 3S3P, 10.9V, 8400mAh. Volume Order (deliver from Taiwan with main unit only)